# ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)

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\*Confidence Intervals

## Limit of Detection

The LOD for ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab) was established using limiting dilutions of a viral sample inactivated. The material (ZeptoMetrix, 0810587CFHI) was supplied at a concentration of  $1.15 \times 10^7$  TCID<sub>50</sub>/mL. The Estimated LOD is 1000 TCID<sub>50</sub>/mL

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W : [www.isotekindo.co.id](http://www.isotekindo.co.id)

## TEST PROCEDURE

Allow the test, specimen, extraction buffer to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed foil pouch and use it within one hour. Best results will be obtained if the assay is performed immediately after opening the foil pouch.
2. Place the Extraction Tube in the workstation. Hold the extraction reagent bottle upside down vertically. Squeeze the bottle and let the solution drop into the extraction tube freely without touching the edge of the tube. Add 8 drops of solution (Approx. 300µL) to the Extraction Tube. See illustration 1.
3. Place the swab specimen in the Extraction Tube. Rotate the swab for approximately 10 seconds while pressing the head against the inside of the tube to release the antigen in the swab. See illustration 2.
4. Remove the swab while squeezing the swab head against the inside of the Extraction Tube as you remove it to expel as much liquid as possible from the swab. Discard the swab in accordance with your biohazard waste disposal protocol. See illustration 3.
5. Fit the dropper tip on top of the extraction tube. Place the test cassette on a clean and level surface. See illustration 4
6. Add 3 drops of the solution (approx.80µL) to the sample well and then start the timer. Read the result at 10 minutes. Do not interpret the result after 20 minutes.

## INTERPRETATION OF RESULTS

(Please refer to the illustrations)

### POSITIVE:

Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T). A positive result indicates that SARS-CoV-2 was detected in the specimen.


**\*NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of SARS-CoV-2 Antigen present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

### NEGATIVE:

One colored line appears in the control line region (C). No line appears in the test line region (T). A negative result indicates that SARS-CoV-2 antigen is not present in the specimen, or is present below the detectable level of the test.

### INVALID:

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Version #: HBT-5354402  
Effective Date: 2020-12

AccuBioTech Co., Ltd.  
Building 10, No. 28 Yu Hua Road, Beijing, China

EC REP Medical Device Safety Service GmbH  
Schiffgraben 41, 30175 Hannover, Germany

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