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

Sinocare

Blood Lipid Reagent Kit (TC/TG/HDL-C/LDL-C)  
TC / TG / HDL-C / LDL-C

13

Blood Lipid Reagent Kit (TC/TG/HDL-C/LDL-C)  
TC / TG / HDL-C / LDL-C

13

**Specifications**

[ Package Size ]  Specification A: 30 tests/kit  
 [ Storage Conditions & Expiry date ]

Product	Expired by	Storage Conditions
Reagent cartridge	2 hours	Room temperature. Before testing, please equilibrate the reagent cartridge to room temperature. (Avoid bright light exposure). Cartridges cannot be returned to the refrigerator once they had been at room temperature.
	12 months	2-8°C.

If applicable

Product	Expired by	Storage Conditions
Control material	Opened and reconstituted	14 days 2-8°C.
		28 days -15~-25°C. The control material can only be frozen once.
	Unopened	12 months 2-8°C.

[ Test environment ]  
 Temperature: 15 °C ~30 °C  
 Relative Humidity: 25%~85%  
 Please perform the test under prescribed conditions, or incorrect results may be obtained.  
 [ Contents ]  Reagent cartridge  Tip  Control material

Blood Lipid Reagent Kit (TC/TG/HDL-C/LDL-C)  
TC / TG / HDL-C / LDL-C

13

**Cautions**  
 1. This product is for in vitro diagnostic use only. Please take protective measures during operation (ex. wearing mask, gloves, etc.). Once contact with eyes or skin, wash immediately with water; go to the hospital if swallow by mistake;  
 2. Please treat used or expired reagent kits as medical waste when disposing them.

CE Sinocare

KEMENKES RI AKL  
 Didistribusikan Oleh:  
 PT Sinocare Healthcare Indonesia  
 Jakarta - Indonesia

(Intended Use) The Blood Lipid Reagent Kit (TC/TG/HDL-C/LDL-C) is intended to quantitatively determine four types of blood lipid: total cholesterol (TC), triglyceride (TG), high-density lipoprotein cholesterol (HDL-C), low-density lipoprotein cholesterol (LDL-C) in human serum. Clinically, it is mainly used as an aid in the diagnosis of hypercholesterolemia, hypertriglyceridemia, coronary heart disease and atherosclerosis.  
 (Applicable Device) Portable Automatic Multi-function Analyzer

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