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

UMONIUM³⁸ INSTRUMENTS

MADE IN BELGIUM

BACTERICIDAL, MYCOBACTERICIDAL, FUNGICIDAL, VIRUCIDAL

Complies with: EN1040:2006, EN1275:2005, EN13727:2003,
EN13624:2004, EN14348:2005, EN14476:2013, EN14561:2006,
EN14562:2006, EN14563:2009

ET0R124 - Ed. 4.2 - 06/2024

BACK

EN Cleaning and disinfection for invasive medical devices before sterilisation. **INSTRUCTIONS:** Cleaning/pre-disinfection: Use at 0.5% (25 mL/5L of water*). Soak. Scrub the devices or use the ultrasounds to eliminate all dirt. Rinse thoroughly with water*. Disinfection: Use at 2.5% (25 mL/L of water*) on clean devices. Make sure to respect a minimum soaking time of 10 minutes. Rinse thoroughly with water*. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

ETOV403 - Ed 1.0 - 05/2024

COMPOSITION UMONIUM³⁸ INSTRUMENTS (SL): N-benzyl-N-dodecyl-N,N-dimethyl-ammonium chloride / N-benzyl-N,N-dimethyl-N-tetradecyl-ammonium chloride. 99g/L.

125 ml e 4,4 oz.fl. **25ml = 2X 12,5ml**

INDONESIA AUTHORIZED REPRESENTATIVE: PT. ISOTEKINDO INTERTAMA,
Jl. Raya Kebayoran Lama No. 309 C, Jakarta 12210, Indonesia
No. Reg. Alat Kesehatan KEMENKES RI AKL xxxxxxxxxxxx

LABORATOIRE HUCKERT'S INTERNATIONAL
20 Avenue Lavoisier - 1300 Wavre - BELGIUM
Tel.: +32 (0)67 89 41 00 - www.huckerts.net

1 L

ETOR298 - Ed. 1.0 - 10/2021

FRONT

UMONIUM³⁸ INSTRUMENTS

CLEANING AND DISINFECTION OF **INVASIVE MEDICAL DEVICES**
PEMBERSIHAN DAN DESINFEKSI **PERANGKAT MEDIS INVASIF**

BACTERICIDAL, MYCOBACTERICIDAL, FUNGICIDAL, VIRUCIDAL

EN1040:2006, EN1275:2005, EN13624:2004, EN13727:2003, EN14348:2005, EN14476:2013, EN14561:2006, EN14562:2006, EN14563:2009

MADE IN BELGIUM
1l e 35.2 fl.oz

UMONIUM³⁸ INSTRUMENTS

BACK

EN INDICATIONS: Cleaning and disinfection for invasive medical devices before sterilisation. **INSTRUCTIONS: Cleaning/pre-disinfection:** Use at 0.5% (25 mL/5L of water*). Soak. Scrub the devices or use the ultrasounds to eliminate all dirt. Rinse thoroughly with water*. **Disinfection:** Use at 2.5% (25 mL/L of water*) on clean devices. Make sure to respect a minimum soaking time of 10 minutes. Rinse thoroughly with water*. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

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COMPOSITION UMONIUM³⁸ INSTRUMENTS (SL): N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride/N-benzyl-N, N-dimethyl-N-tetradecyl-ammonium chloride. 99g/L.

INDONESIA AUTHORIZED REPRESENTATIVE: PT. ISOTEKINDO INTERTAMA.
JL. Raya Kebayoran Lama No. 309 C, Jakarta 12210, Indonesia
MEDICAL DEVICE CLASS B - No. Reg. Alat Kesehatan KEMENKES RI AKL XXXXXXXXXXXX

LABORATOIRE HUCKERT'S INTERNATIONAL
20 Avenue Lavoisier - 1300 Wavre - BELGIUM
Tel.: +32 (0)67 89 41 00 - www.huckerts.net


ETOV293 - Ed. 1.0 - 10/2021

5 L

FRONT

ET0R299 - Ed. 1.0 - 10/2021

UMONIUM³⁸ INSTRUMENTS

CLEANING AND DISINFECTION OF INVASIVE MEDICAL DEVICES
PEMBERSIHAN DAN DESINFEKSI PERANGKAT MEDIS INVASIF

BACTERICIDAL, MYCOBACTERICIDAL, FUNGICIDAL, VIRUCIDAL

EN 1040:2006, EN 1275:2005, EN 13624:2004, EN 13727:2003, EN 14348:2005, EN 14476:2013, EN 14561:2006, EN 14563:2009

MADE IN BELGIUM
5l e 176 fl.oz

UMONIUM³⁸ INSTRUMENTS

BACK

EN INDICATIONS: Cleaning and disinfection for invasive medical devices before sterilisation. **INSTRUCTIONS: Cleaning/pre-disinfection:** Use at 0.5% (25 mL/5L of water*). Soak. Scrub the devices or use the ultrasounds to eliminate all dirt. Rinse thoroughly with water*. **Disinfection:** Use at 2.5% (25 mL/L of water*) on clean devices. Make sure to respect a minimum soaking time of 10 minutes. Rinse thoroughly with water*. **Clean drinking water at room temperature.* **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

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COMPOSITION UMONIUM³⁸ INSTRUMENTS (SL): N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride / N-benzyl-N, N-dimethyl-N-tetradecyl-ammonium chloride. 99g/L.

INDONESIA AUTHORIZED REPRESENTATIVE: PT. ISOTEKINDO INTERTAMA.
Jl. Raya Kebayoran Lama No. 309 C, Jakarta 12210, Indonesia
MEDICAL DEVICE CLASS B - No. Reg. Alat Kesehatan KEMENKES RI AKL XXXXXXXXXXXX

LABORATOIRE HUCKERT'S INTERNATIONAL
20 Avenue Lavoisier - 1300 Wavre - BELGIUM
Tel.: +32 (0)67 89 41 00 - www.huckerts.net

ET0V293 - Ed. 1.0 - 10/2021