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
PH
NEUTRAL

MEDICAL
DEVICE
Class IIb

EN Cleaning and pre-disinfection for invasive medical devices before sterilisation. Prepare for disinfection. **INSTRUCTIONS: Cleaning/Pre-disinfection:** Use at 0.5% (25 mL/5L of water *). Soak then scrub/use the ultrasounds to eliminate all dirt. Rinse thoroughly with clean water before disinfection. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

BACK

ETOV402 - Ed 1.0 - 05/2024

COMPOSITION UMONIUM³⁸ **STERILY** (SL): N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride/N-benzyl-N,N-dimethyl-N-tetradecyl-ammonium chloride. 22 g/L

125 ml e 4,4 oz.fl. **25ml = 2X 12,5ml**

INDONESIA AUTHORIZED REPRESENTATIVE: PT. ISOTEKINDO INTERTAMA.
Jl. Raya Kebayoran Lama No. 309 C, Jakarta 12210, Indonesia
No. Reg. Alat Kesehatan KEMENKES RI AKL XXXXXXXXXXX


LABORATOIRE HUCKERT'S INTERNATIONAL
20 Avenue Lavoisier - 1300 Wavre - BELGIUM
Tel: +32 (0)67 89 41 00 - www.huckerts.net

1 L

ETOR302 - Ed 1.0 - 10/2021

FRONT

UMONIUM³⁸ STERILY

CLEANING AND PRE-DISINFECTION FOR INVASIVE MEDICAL DEVICES
PEMBERSIHAN DAN PRADISINFEKSI UNTUK PERANGKAT MEDIS INVASIF

MADE IN BELGIUM
1l e 35.2 fl.oz

UMONIUM³⁸ STERILY

BACK

EN INDICATIONS: Cleaning and pre-disinfection for invasive medical devices before sterilisation. Prepare for disinfection. **INSTRUCTIONS: Cleaning/Pre-disinfection:** Use at 0.5% (25 mL/5L of water *). Soak then scrub/use the ultrasounds to eliminate all dirt. Rinse thoroughly with clean water before disinfection. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

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COMPOSITION UMONIUM³⁸ STERILY (SL): N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride/N-benzyl-N, N-dimethyl-N-tetradecyl-ammonium chloride. 22 g/L

INDONESIA AUTHORIZED REPRESENTATIVE: PT. ISOTEKINDO INTERTAMA.
Jl. Raya Kebayoran Lama No. 309 C, Jakarta 12210, Indonesia
MEDICAL DEVICE CLASS B - No. Reg. Alat Kesehatan KEMENKES RI AKL XXXXXXXXXXXX

LABORATOIRE HUCKERT'S INTERNATIONAL
20 Avenue Lavoisier - 1300 Wavre - BELGIUM
Tel.: +32 (0)67 89 41 00 - www.huckerts.net

ET0V295 - Ed. 1.0 - 10/2021

5 L

FRONT

ETOR303 - Ed 1.0 - 10/2021

UMONIUM³⁸ STERILY

CLEANING AND PRE-DISINFECTION FOR **INVASIVE MEDICAL DEVICES**
PEMBERSIHAN DAN PRADISINFEKSI UNTUK **PERANGKAT MEDIS INVASIF**

MADE IN BELGIUM
5l e 176 fl.oz

UMONIUM³⁸ STERILY

BACK

EN INDICATIONS: Cleaning and pre-disinfection for invasive medical devices before sterilisation. Prepare for disinfection. **INSTRUCTIONS: Cleaning/Pre-disinfection:** Use at 0.5% (25 mL/5L of water *). Soak then scrub/use the ultrasounds to eliminate all dirt. Rinse thoroughly with clean water before disinfection. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

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