

Seasight ESR Tube



Tabung dan Pipet Ukur


- Material : USP VI Polypropylene
- Skala baca : 0-150mm
- Cairan antikoagulan : 0.2 ml (sodium citrate 3,8%)
- Waktu baca : 60 Menit
- Metode : Westergren
- ED : 2 Tahun

Rak

- Material : Polypropylene
- Jumlah channel : 10 ruang



Seasight ESR Tube 170nm



Tabung dan Pipet Ukur

- Material : USP VI Polypropylene
- Skala baca : 0-170mm
- Cairan antikoagulan : 1,6 ml (Na citrate 3,8%)
- Waktu baca : 60 Menit
- Metode : Westergren
- ED : 2 Tahun

Rak

- Material : Polypropylene
- Jumlah channel : 10 ruang

Seasight Fast Detector



Tabung dan Tabung Ukur

- Material : USP VI Polypropylene
- Skala baca : 0-140mm
- Cairan antikoagulan : 1,28 ml (Na citrate 3,8%)
- Waktu baca : 30 Menit
- Metode : Fast Detector westergren
- ED : 2 Tahun

Rak

- Material : Polypropylene
- Jumlah channel : 10 ruang

